


Advazorb

Hydrophilic Foam Dressing Range



Gentle on skin

Distributed by



www.mediusa.com

Adv  **ancis**
Medical

www.advancismedical.com

Advazorb®

Hydrophilic foam dressing range

Manufactured in the UK, the Advancis Medical® range of Advazorb® Hydrophilic foam dressings, offer the combination of a high quality performance product at a competitive price.

Advazorb® is a comprehensive range of patient friendly, absorbent foam dressings presented in non-adhesive and atraumatic silicone adhesive format. The range is designed to overcome the complex challenges of managing exudate whilst protecting 'at risk' fragile skin, a combination of excellent fluid handling capabilities of Advazorb® with a layer of atraumatic soft silicone adhesive.

Advazorb®, Advazorb Silfix® and Advazorb Border® are available in both a regular thickness for moderate to high exudate and 'Life' versions for low to moderately exuding wounds.



Advazorb Border®

Hydrophilic foam dressing with soft silicone wound contact layer and border

The silicone wound contact layer has large pores which enable the passage of exudate into the absorbent foam whilst protecting the fragile wound bed.

This combination ensures the dressings comfortably stay in place whilst minimizing the potential for maceration.

- The same high quality foam and low friction film as Advazorb® and Advazorb Silfix®, plus a border to eliminate the need for secondary fixation.
- The soft silicone contact layer extends across the entire surface area of the dressing, not just the border, providing secure but gentle adhesion.
- Suitable for hard to dress wounds as a secondary fixation isn't required.
- The large open pores enable the passage of exudate through to the foam whilst providing optimal adherence to the skin.
- Minimizes the pain and trauma associated with dressing change or removal.



Advazorb Silfix®

Hydrophilic foam dressing with soft silicone wound contact layer

Combination of the excellent fluid handling capabilities of Advazorb® with a layer of atraumatic soft silicone adhesive.

The soft silicone ensures that the dressing does not stick to the wet surface of a wound bed or cause trauma to delicate new tissue upon dressing removal. This is particularly suited to those patients with friable or delicate skin.

- The same high quality foam and low friction film as Advazorb® with an atraumatic soft silicone wound contact layer across the surface area of the dressing.
- The wound contact layer enables the dressing to be placed onto the patient whilst a secondary dressing or bandage is applied for retention.
- Provides gentle adherence to intact skin whilst the hydrophobic nature of silicone ensures that it won't adhere to a wet wound bed.
- Minimizes the pain and trauma associated with dressing change or removal.
- Prevents disruption of newly formed granulation tissue.
- Ideal choice for wounds with fragile skin that also require a degree of compression.

Product	Size	Dressings per box	HCPCS code
Advazorb Border®	3in x 3in	10	A6212
	4in x 4in	10	A6212
	5in x 5in	10	A6213
	6in x 6in	10	A6213
	4in x 8in	10	A6213
	8in x 8in	10	A6214

Product	Size	Dressings per box	HCPCS code
Advazorb Silfix®	3in x 3in	10	A6209
	4in x 4in	10	A6209
	5in x 5in	10	A6210
	6in x 6in	10	A6210
	4in x 8in	10	A6210
	8in x 8in	10	A6211



medi USA
6481 Franz Warner Pkwy
27377 Whitsett
United States of America
Tel: 800-633-6334 | E-mail: info@mediusa.com

MAR437/IV2



www.advancismedical.com